

**zygo**


---

**Basic Operator Procedure Training**  
**NewView 6300**




# NewView 6200/6300

- System Overview



# NewView Theory of Operation

- Measurement Technology: Scanning White-light Interferometry (SWLI)
  - Interferometric objective mounted in a piezo scanning device that moves vertically
  - Camera detects interferogram intensities
  - Computer stores only good (modulating) data as 3D interferogram
  - Frequency domain analysis (FDA) determines heights of each pixel to  $<0.1\text{nm}$  resolution and  $<0.3\text{ nm}$  repeatability




# Motion Controller/Joystick

Set Z STOP  
To avoid crashing  
objective into  
sample

Change Z AXIS  
SPEED to speed up  
or slow down travel  
for focusing

Use AXIS  
SELECT to select  
X Y for moving  
stage or R P for  
adjusting tip/tilt

FOCUS JOYSTICK  
Move up to raise z  
away from stage;  
Move down to lower  
toward stage



Press MOTION  
STOP to halt  
motion  
in case of  
emergencies


Move STAGE JOYSTICK left/right to  
move in X or R and up/down to move  
in Y or P

# Operator Process Flow




## Load Application

- Open the MetroPro Shortcut on the Desktop
- Click on the desired application to open



# Load Application

- If application not on desktop, right-click to “Load Application”
- Select and load application from list



# Typical MetroPro Layout

- Micro.app: generic application for measuring surface structure and roughness of samples

The screenshot shows the Zygo MetroPro software interface with the following components labeled:


- Objective Button**: Points to the '5X SLWD' button in the top left.
- Buttons**: Points to the 'MEASURE' menu containing 'Analyze', 'Mask Data', 'Save Data', 'Load Data', 'Calibrate', and 'Reset'.
- Control Windows**: Points to the 'Measure Cntrl', 'Analyze Cntrl', 'Stage Controls', 'Focus Controls', and 'Sequence Controls' windows.
- Filled Plot**: Points to the 'Surface Map' window showing a color-coded surface topography.
- Numerical Results**: Points to the 'Measure Attributes' window at the bottom, which displays:
 

|                          |                                 |
|--------------------------|---------------------------------|
| Fri Apr 28 15:38:43 2006 | Objective: 20X Mirau            |
| P/N:                     | Subtract Sys Err: Off           |
| S/N:                     | Camera Res: 0.552 $\mu\text{m}$ |
- Surface Profile**: Points to the 'Surface Profile' window showing a line graph of Height ( $\mu\text{m}$ ) vs Distance (mm).
- 3D Model**: Points to the '3D Model' window showing a 3D surface topography plot.
- Solid Plot**: Points to the 'Solid Plot' window showing a grayscale surface texture image.
- Measure Attributes**: Points to the 'Measure Attributes' window at the bottom right.



## Select Objective

- Click the Objective Button to select the objective
- For automated turret, ensure sample stage is clear first
- Choose objective based on resolution and desired field of view
  - Tradeoff: Lower magnification yields a wider field of view but provides coarser lateral resolution



## Load Sample and Set Z-Stop


- Place test sample on stage under objective
- Using joystick, lower head to slightly less than the objective working distance
- Press “SET” Z-stop
  - Important step to avoid crashing objective into sample

| Objective             | 1X  | 2.5X | 5X  | 10X | 20X | 50X | 100X | SLWD |
|-----------------------|-----|------|-----|-----|-----|-----|------|------|
| Working Distance (mm) | 8.3 | 10.3 | 9.3 | 7.4 | 4.7 | 3.4 | 0.8  | 40   |

## Focus Sample

- Adjust z using the focus joystick until sample is focused
- Focus is found when fringes appear on Live Display
  - Fringes are the light and dark bands produced by the interference of light
- *Hints:* if having trouble focusing, press “F5” to set the light level automatically for viewing, or try using the field stop

Example: Part focused with fringes




## Focus Sample: Field Stop Focus Aid

- The field stop on the side of the NewView can be used as a focus aid
- When the sample gets close to focus, a second image of the iris will appear blurred; when they exactly overlap, fringes should appear
- Most effective for high magnifications and smooth samples

System Out of Focus




System In Focus




## Null Fringes

- Nulling is the process of minimizing the number of fringes
- Null sample by adjusting tip/tilt or “R P” using the joystick
- As the tip/tilt is adjusted, it may be necessary to make fine z-adjustments to keep the fringes viewable on the live display monitor

Fringes due to Tilt  
NOT NULLED




Minimum Fringes  
NULLED




# Null Fringes

- *Note:* The null fringe location will look different depending on the part
  - in most cases, think of “spreading out” the fringes
  - for a spherical part, center the bulls eye


## Example of Nulled Parts



**Smooth Flat Part**  
Adjust for high contrast and the least number of fringes.



**Rough Flat Part**  
The fringes are in smaller isolated areas. Center the fringes and adjust focus between the high and low fringes.




**Spherical Part**  
Adjust the stage and focus to center the circular fringe pattern.

# Set Light Level

- Set **automatically** by pressing “F5”
  - Must be centered on the brightest fringe
- Set **manually** by pressing “F4”
  - Use numeric keypad to set peak intensity to approx 90 – 99 %
  - Make sure there is no saturation (red)

| Key  | Function   |
|--|--|
| <input type="button" value="/"/> <input type="button" value="*"/>                                      | Coarse down (/) and up (*).                          |
| <input type="button" value="-"/> <input type="button" value="+"/>                                      | Fine down (-) and up (+).                            |
| <input type="button" value="1"/> <input type="button" value="2"/> ... <input type="button" value="9"/> | Adjust level in 10% steps, 1 = 10%, 3 = 30%, 9 = 90% |
| <input type="button" value="Tab"/>   | Toggle between last two light settings.              |
| <input type="button" value="Esc"/>   | Cancel current action and abort measurement.         |



# Set Measurement Controls



The screenshot shows the 'Measurement Controls' window in the Zygo software. The window is divided into several sections for configuring measurement parameters. Annotations with arrows point to specific fields:

- User entry for part identification:** Points to the 'Part Number' and 'Serial Number' fields.
- Adjust to acquire more data points:** Points to the 'Min Mod (°): 15' and 'Min Area Size: 7' fields.
- Using averaging to improve accuracy:** Points to the 'Number of Averages: 0' field.
- Use to subtract system errors:** Points to the 'Subtract Sys Err: Off' and 'Sys Err File: SysErr.dat' fields.
- For use with extended scan option:** Points to the 'Extended Scan Length: 1000 um' field.

The interface includes the following controls:

- Comment:** A text area for notes.
- Part Number:** A text field for the part's ID.
- Serial Number:** A text field for the part's serial number.
- Min Mod (°): 15** and **Min Area Size: 7**: Numerical input fields for scan parameters.
- Instrument:** A dropdown menu currently set to 'None'.
- Acquisition Mode:** A dropdown menu set to 'Scan'.
- Camera Mode:** A text field showing '640x480 22Hz'.
- Image Zoom:** A dropdown menu set to '1X'.
- Remove Fringes:** A dropdown menu set to 'On'.
- Subtract Sys Err:** A dropdown menu set to 'Off'.
- Sys Err File:** A text field containing 'SysErr.dat'.
- Number of Averages:** A numerical input field set to '0'.
- Phase Controls:**
  - AGC:** A dropdown menu set to 'On'.
  - Phase Res:** A dropdown menu set to 'High'.
  - Connection Order:** A dropdown menu set to 'Location'.
  - Discon Action:** A dropdown menu set to 'Filter'.
- Scan Controls:**
  - Scan Length:** A dropdown menu set to 'None'.
  - Extended Scan Length:** A text field set to '1000 um'.
  - FDA Res:** A dropdown menu set to 'None'.



## Measurement Controls


- Image Zoom to increase or decrease field of view, if available
  - 1X standard; 0.5X, 0.75X, 1.5X and 2.0X optional
- FDA Res control sets how the software processes the data collected
  - “High 2G” for smooth surfaces
  - “Normal” for rough surfaces, typically > 75 nm Ra
  - “Low” required for extended scans
- Camera Mode
  - Selects effective camera size for collecting data
  - More pixels resolve smaller details but result in increased processing time

# Measurement Controls


## • Scan Length

- Select length of measurement scan
- Ranges from 5  $\mu\text{m}$  to 15 mm
- Longer scan = Longer acquisition time
- Bipolar Scan: from initial position, objective moves down half the scan length and then scans upward
- Extended Scan scans upward only

### Bipolar Scan




### Extended Scan



## Measurement Controls

- Min Mod (%)
  - Specifies minimum modulation or intensity range for a valid data point
  - Setting can range from 0 to 100 %



# Analyze Controls

Removes "form"  
from the data

Removes  
edge effects

Fill holes in  
the data

Filtering controls  
are used to select  
and analyze the  
various spatial  
frequency  
components that  
make up the  
test part.

zygo Analyze Controls

|                         |                           |
|-------------------------|---------------------------|
| Remove: Plane           | Min Area Size: 0          |
| Sphere Radius: 0 nm     | Remove Spikes: Off        |
| Trim: 0                 | Spike Height (xRMS): 7.50 |
| Trim Mode: All          |                           |
| Data Fill: Off          |                           |
| Data Fill Max: 25       |                           |
| Filter: Off             |                           |
| Filter Type: Average    |                           |
| Filter Window Size: 3   |                           |
| Filter Trim: Off        |                           |
| Filter Low Wavelen: mm  | Filter Low Wavelen: nm    |
| Filter High Wavelen: mm | Filter High Wavelen: nm   |
| Filter Low Freq: 1/mm   | Filter Low Freq: 1/mm     |
| Filter High Freq: 1/mm  | Filter High Freq: 1/mm    |

- *Note:* Changes to these controls can be made after measurement  
Press "F2" or **Analyze** button to re-analyze after changes


## Analyze Controls

- Remove specifies the surface to remove to minimize form
  - As a general rule of thumb, remove a plane for flatness or a cylinder for roughness
- Turn Data Fill On to fill missing data points; The maximum number of pixels that will be filled is based on the Data Fill Max control
- If Remove Spikes is On, a pixel will be removed if its height is greater than the surrounding pixel heights by the Spike Height value

# Analyze Controls

## • Filtering



- Low pass, high pass, band pass or band reject filters are available in the Filter control (Off by default)
  - Low pass highlights waviness or form; high pass highlights roughness
- Use Filter Type to choose an average, median, 2 sigma, FFT or Gaussian type filter
- For FFT fixed, enter cutoff values in the high and low wavelength (or frequency) controls



## Measure

- Press “F1” or click  button to measure
- System will scan then display results
  - Do not touch vibration isolation table or sample stage until measurement is complete



# Numerical Results

| Result        | Description   |  |
|---------------|---|--|
| <b>PV</b>     | The height between the lowest and the highest point on the test part surface.                             |  A line graph showing a jagged surface profile. Two horizontal lines are drawn: one at the highest peak and one at the lowest valley. A vertical double-headed arrow between these two lines is labeled 'PV'.                                   |
| <b>rms</b>    | The root-mean-square deviation of all points from a plane fit to the test part surface.                   |  A line graph showing a jagged surface profile. A solid horizontal line represents the plane fit, and a dashed horizontal line represents the mean. A vertical double-headed arrow between the solid line and the dashed line is labeled 'rms'. |
| <b>Ra</b>     | The average roughness, or the average deviation, of all points from a plane fit to the test part surface. |  |
| <b>Size X</b> | The dimension of the data set in the X-axis (on-screen horizontal).                                       |  |
| <b>Size Y</b> | The dimension of the data set in the Y-axis (on-screen vertical).   |  |



## Save Data: Store screenshot


- Options to Save Data: Save a screenshot, save the raw data or save process stats
- To acquire a screenshot from MetroPro
  1. Click on the Zygo button on the upper left corner of the application
  2. Choose File, .bmp, Color from the Print Panel then click Print
  3. Save the file with a .bmp extension



# Save Results: Save Process Stats

- Click the process icon to open the process window
- To save a series of data as a txt or csv file, press the zygo button and choose to print to a file

The Process Window



Store- adds one set of data  
 Undo- removes last data set  
 Clear- erases ALL data sets

Window Size sets the number of rows displayed.

Each row shows the results for one measurement. Rows fill in as measurements are made when Auto Store is On.

|        | PV<br>um | rms<br>um | Ra<br>um |
|--------|----------|-----------|----------|
| 1      | 10.089   | 0.915     | 0.730    |
| 2      | 9.171    | 0.898     | 0.717    |
| 3      | 6.826    | 0.126     | 0.070    |
| 4      | 10.089   | 0.915     | 0.730    |
| 5      |          |           |          |
| 6      |          |           |          |
| 7      |          |           |          |
| 8      |          |           |          |
| 9      |          |           |          |
| 10     |          |           |          |
| Range  | 3.263    | 0.790     | 0.660    |
| Mean   | 9.044    | 0.714     | 0.562    |
| StdDev | 1.540    | 0.392     | 0.328    |
| NValid | 4        | 4         | 4        |

Range- difference between largest and smallest values.  
 Mean- arithmetical average of all values.  
 StdDev- standard deviation  
 NValid- number of valid values

## Save Results: Data

- Press `Save Data` button
- In the file handler, enter a name for the file ending with “.dat”
- Raw data is saved; Can be post-analyzed

## Troubleshooting

- What to do if MetroPro returns “No Valid Data” error or if there is too much data dropout
  - Check Focus: Do you see fringes and are they nulled? (Or for an extended scan, are you positioned below focus?) 3 or fewer fringes is a good rule of thumb
  - Check Scan Length: Is your scan long enough? Do you see all of the fringes go by?
  - Check Light Level (F4): Is light level in green zone?
  - Check Min Mod %: Does it need to be lowered?

## Resources

- NewView 6200/6300 Operating Manual, OMP-0503
- NewView Microscope Application Booklet, OMP-0360
- MetroPro Reference Guide, OMP-0347
- Available for download online:

<http://www.zygo.com/?sup=/resource/manuals.cgi>

User ID: manuals

Password: d5g2r8

*Note:* User ID and password periodically changing; Contact Zygo if unable to access